




186D FISHLINE INKER



The 186D Inker mechanism is used to mark individual die during the testing process. The finger adjustable controls are all conveniently located on top of the Inker for easy access. Inkers can be cleaned frequently as ink reservoirs are easily removed without disturbing other parts. Varied fishline sizes are available to create different ink dot sizes.	APPLICATION	The 186D Inker mechanism is used to mark individual die during the testing process.				
	PHYSICAL DIMENSIONS	Width	4"			
		Length	4"			
Height		3.5"				
AVAILABLE ACCESSORIES	Dot/Pack Size (mils)	Inker	Head Assembly	Reservoir	Monofilament	
	Less reservoir & filament	17181-000	11780	-	-	
	5 (soft)	17181-005	11780-005	11705-004	14101	
	5 (hard)	17181-105	11780-105	11705-005	11708	
	10	17181-010	117080-010	11710-003	11707-006	
	15	17181-015	11780-015	11715-003	11712-006	
	20	17181-020	11780-020	11720-003	11717-006	
	25	17181-025	11780-025	11725-003	11722-006	
	30	17181-030	11780-030	11730-003	11727-006	
	Inker Manipulator, Bracket and Magnetic plate	17181-050	-	-	-	

FOR ORDERING INFORMATION CONTACT SALES (408) 934-9875